



Compatibility Chart Mac OS X

OS version	AM211/AM310/AM311	AM311S(T)	AM313X series	413/423 series	3613/3713 series	2011/3013/4013/4023 series	7013/7023 series
Mac OS 9.1 and up	yes	yes	yes	yes	no	no	no
Mac OS X 10.0	yes	yes	yes	yes	no	no	no
Mac OS X 10.1	yes	yes	yes	yes	no	no	no
Mac OS X 10.2	yes	yes	yes	yes	no	no	no
Mac OS X 10.3	yes	yes	yes	yes	no	no	no
Mac OS X 10.4	yes	yes	yes	yes	no	yes	yes
Mac OS X 10.5 (Leopard)	no	yes	yes	yes	no	yes	yes
Mac OS X 10.6 (Snow Leopard)	no	yes	yes*	yes	no	yes	yes
Mac OS X 10.7 (Lion)	no	yes	yes*	yes	no	yes	yes

-  DinoXcope/DinoScope
-  DinoScope Only
-  DinoXcope Only

* Only on systems that COME with that OS, not for upgrades